




BETTER HISTOLOGY SOLUTION
BETTER LIFE FROM HEALTH



ASD 220 SLIDE DRYER

Technical Specifications

- Temperature control accuracy: $\pm 1^{\circ}\text{C}$
- Temperature range: ambient temperature $\sim 90^{\circ}\text{C}$
- Capacity : can load 60 slides simultaneously
- Rated voltage:
220V $\pm 10\%$, 50/60 Hz / 110V $\pm 10\%$, 50/60 Hz
- Nominal power: 330W
- Outside dimensions(mm): 330(W) \times 305(D) \times 85(H)
- Weight: 2.75kg



AMOS SCIENTIFIC PTY.LTD
Address: 37/632 CLAYTON ROAD ,
CLAYTON SOUTH VIC 3169 , AUSTRALIA
Tel: +61-3-95518850
Fax: +61-3-95516885
Website: www.amos-scientific.com
Email : info@amos-scientific.com