



BETTER HISTOLOGY SOLUTION
BETTER LIFE FROM HEALTH



TEC2601 WATER BATH

Technical data of TEC2601 Water Bath

- Water trough dimensions: 230 × 180 × 56 mm
- Temperature control accuracy: $\pm 1^{\circ}\text{C}$
- Temperature range: ambient temperature $\sim 75^{\circ}\text{C}$
- Rated voltage: 220V $\pm 10\%$ 50/60 Hz or 110V $\pm 10\%$ 50/60 Hz
- Nominal power: 350W
- Outside dimensions (W × D × H): 328 × 300 × 85 mm
- Weight: about 2.4kg
- Capacity : 2.3 L

