

## Cryo Console

- The paraffin block can be cold down rapidly by the stable and environmental friendly compressor
- The temperature range is between -15 to 15° C
- 70 cassettes can be loaded on this large work surface at the same time
- User-friendly operation interface



## Cryo Console

Voltage/Frequency	220 V ± 10% , 50 HZ / 110 V ± 10% , 60 HZ
Nominal power	330 W
Working temperature	15°C~-15°C    Increments 1°C

### Dimensions

Width	+	390mm
Depth		560mm
Height		390mm
Weight		about 24 KG
Capacity		70 cassettes

