



EC500-2: DISPENSING CONSOLE

Dimensions (w x d x h)

Console	345 x 641 x 405 mm
Heated area	340 x 253 mm

Weight 20,5 Kg

Adjustable temperature range

Paraffin reservoir	40-70 °C
Heated area/ forceps wells	40-70 °C
Electrically heated tweezers	60-75 °C

Paraffin reservoir capacity 5 l.

Cold spot temperature 5 °C

Accessories Magnifier

Footswitch

Electrically heated tweezers, 1, 2 and 4 mm tip



EC500-3: THERMAL CONSOLE

Dimensions (w x d x h)

Console	340 x 623 x 395 mm
Cassette storage	269 x 510 x 88 mm
Removable trays	262 x 175 x 55 mm
Mold storage	275 x 170 x 180 mm

Weight 16 kg

Adjustable temperature range (via Dispensing Console)

Cassette storage	40-70 °C
Mold storage	40-70 °C



EC500-1: CRYO CONSOLE

Dimensions (w x d x h)

Console	335 x 623 x 395 mm
Cold plate	375 x 307 mm

Weight 24,5 kg

Adjustable temperature range de 0 a -12 °C

Especialidades Médicas Myr S.L

ISO 9001 / 13485 certified company

C / Vidriers, 19 - 21, Pol. Ind. La Cometa 43700

El Vendrell, Spain

Phone. +34 977 66 8020

Fax. +34 977 66 8030

Mail. esp.medicas@myr.com.es

www.myr.com.es

Local Distributor. Contact information.