

# Laser Cassette Printer

## PrintSmart LCP-280



High Resolution




Permanent Storage



Orderliness Collection



Fast Printing



Batch Collection



Automatic Cleaning




### Specifications

<b>Capacity of loading</b>	6 Hoppers for 600 Cassettes	<b>Data interface</b>	USB
<b>Capacity of collection</b>	288 cassettes	<b>Working temp Range</b>	5°C-40°C
<b>Print Speed</b>	3 seconds/cassette	<b>Power Draw</b>	120VA
<b>Print Resolution</b>	2500dpi	<b>Voltage</b>	AC100-240v/50-60hz
<b>Type of cassette</b>	Regular cassette with 45°angles	<b>Dimension</b>	463 x 531x496mm(WxDxH)
<b>Data Acquisition</b>	HIS or LIS	<b>Weight</b>	46kg(net) / 67kg(gross)

# Choosing us

- ✓ Laser Printing on Cassette ,No ink cartridge or ribbon is required which greatly reduces the operation cost.
- ✓ 6 hoppers which allow cassettes in different color and set groups.600 cassettes max for one loading.
- ✓ Durability , Reliable and Quick.
- ✓ Tray collection with large collection quantity,Place 9 trays at a time,and collect up to 288 cassettes totally.
- ✓ High clear print resolution in 0.5mm characters and Bar code is always excellent with good legibility
- ✓ Translucent viewing window allows to observe the printing process at any time, and easy to service.
- ✓ Allow to print characters: any language,texts,symbol,data, logo. barcode etc.,and inserted image.
- ✓ Unlimited the number of print fields ,content, and adjust the print grayscale.
- ✓ Printing Style: Single,Batch, Bar code scanning, Compatible to LIS and HIS online.
- ✓ Friendly operation interface is simple and easy to understand and edit.
- ✓ Power off protection function, Safe and Flexibility.
- ✓ Permanent Storage, Resistant to hard objects and any reagent.
- ✓ The automatic pipeline cleaning function keeps the air fresh and protect the health of operators.
- ✓ Language: English or others as customer demand.

 [info@tria-link.com](mailto:info@tria-link.com)

 **Tel(whatsapp):**  
**+86-013382819886**

 **CHANGZHOU TRIALINK CO., LTD**  
No1408 Building 3,Manhattan Plaza,  
Xinbei District Changzhou City, Jiangsu, China

**TriaLink Co.,Ltd(UK)**  
Unit G25 Waterfront Studios, 1 Dock Road,  
London. E16 IAH,UK

