



PrintSmart LSPlus-X



PrintSmart LSPlus-S



PrintSmart LSPlus-PAD



PrintSmart LASER SLIDE PRINTER

features

microscope slides laser labeler

It prints directly indelible identification information on the slide to create permanent.

Ink cartridges or ribbons dont needed, the operation cost is greatly reduced when compared with the traditional printers.

It can be used for printing characters : texts,symbol,rafphics, logos, barcodes for reliable specimen tracking.

Durability and Reliable

- Clear print quality and permanent Labels for histology process well done and many years in storage
- Reliable Identification to enable the best tracking capacity in your laboratory.

Efficient issues:

- A staggering amount of information on Slides printing to match your laboratory' s need.
- High Clear print resolution is always excellent with good legibility.
- Reduced Human Errors and flexibility.
- quick print speed provides effective and efficient work

keys:

- Loading Capacity: single slot holds up to 100 slides.
- Touch Screen with 8 inches, easy and friendly operation.
- Compatible to LIS and HIS
- Language: English or others as customer demand.

Technical parameter

Capacity of loading :	Single Slot for 100 slides	Data Acquisition:	barcode reader
Capacity of Collection Tray	40slides	Data interface:	WIFI R45
Print Speed:	4 seconds/slide More than 15slides/Minute	Working Temp. Range	5°C-40°C
Print Resolution:	2500dpi	Power Draw	120VA
Type of Slide:	26x76mm, Standard	Voltage:	AC100-240v/50-60hz
Weight:	20kg(net)/45kg(gross)	Dimensions:	220 x 488 x357mm(WxDxH)



TRIA-LINK

Changzhou Trialink Co.ltd

Add:1408, Building 3, Manhattan Plaza, Xinbei District, Changzhou City, Jiangsu, China

Website: <http://www.tria-link.com> Email : info@tria-link.com

Trialink Co.ltd

Add: Unit G25 Waterfront Studios, 1 Dock Road, London,UK